# HLA B5701 STRIP

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
## Detection of HLA B5701 alleles

Abacavir is used in the treatment of HIV-1 infection since 1999. The WHO recommends abacavir as a second-line treatment, owing to the risk of hypersensitivity associated with its use, which affects 5-8% of Caucasian patients.

Several studies have demonstrated the existence of a strong predictive correlation, in Caucasian populations and groups of Hispanic ethnicity, between hypersensitivity to abacavir and the presence of the HLA-B\*5701 allele. This correlation is strong enough to enable us to predict the risk of hypersensitivity to abacavir and classify individuals as low (<1%) or high (>70%) risk, based on the absence or presence of the HLA-B\*5701 allele.

Several polymorphisms has been identified. These ones explain nearly 15% of the variability in HIV-1 viral load in individuals during the asymptomatic phase of the disease were identified. One of them, located in the HCP5 gene (HLA complex P5), shows a perfect linkage disequilibrium with HLA-B\*5701 ( $r^2 = 1$ ).

The useful role of this SNP in predicting hypersensitivity to abacavir has been demonstrated in European populations (100% detection sensitivity for HLA-B\*5701, specificity of 99.4%, positive predictive value of 94.2% and negative predictive value of 100%).

## Intended use

HLA B5701 Strip is a test designed for the detection of HLA B5701 alleles associated with Abacavir hypersensitivity reaction.

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
\* OpeGen kits include all the necessary reagents even Taq polymerase


## RESULTS

### Resultados

A B\*5701 positive sample must always be positive for the bands S1 and S2 as well as for the bands associated with the two control genes (Control A and Control B). All results in which the bands corresponding to the two amplification controls fail to appear must be regarded as invalid.

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- 1) Saag M et al. "High sensitivity of human leukocyte antigen-B\*5701 as a marker for immunologically confirmed abacavir hypersensitivity in white and black patients". *Clinical Infectious Diseases* (2008); 46: 1111-1118.
- 2) Mallal M.B. et al. "HLA-B\*5701 screening for hypersensitivity to Abacavir". *The New England Journal of Medicine* (2008); 358 (6): 568-579
- 3) Colombo S. et al. "The HCP5 single-nucleotide polymorphism: a simple screening tool for prediction of hypersensitivity to reaction to abacavir". *JID* (2008); 198: 864-867.